

Hi-Tech Green 2X qPCR Universal Mix

Cat. No. GQM-M-002-100

Size: 1ml (100 Reactions)

For shipping can be store at 2-4°C for two weeks,
Store at -20°C for long-term under dark condition

MOLEQULE-ON[®]

Kit Components	Quantity
Hi-Tech Green 2X qPCR Universal Mix	1 ml


Description

Hi-Tech Green 2X qPCR Universal Master Mix is a 2X concentrated, ready for use Master Mix reaction enhanced for dye-based quantitative PCR (qPCR) and compatible with the majority of commercially available real-time PCR systems (ROX-independent and ROX-dependent). It contains NanoTaq hot-start DNA polymerase, dNTPs, MgCl₂, Hi-Tech Green I dye, enhancers, stabilizers and essentials for a success PCR reaction.

Features & Applications

Hi-Tech Green qPCR Mix is ideally suited for:

- Enhanced nano complex-mediated hot start enzyme DNA polymerase.
- Maximizes product yield with nano technology.
- Nonspecific bands are eliminated during reaction.
- Compatible with the majority of qPCR systems.
- Gene expression analysis (absolute and relative)
- Detection of low copy genes
- Quantification of viral loads or NGS libraries



		Ct
A	Hi-Tech Green 2X qPCR Universal Mix	21.29
B	Competitor 1	24.70
C	Competitor 2	23.75

Performance of Hi-Tech Green 2X qPCR Universal Mix

A 500 bp human genomic DNA target was used to compare the two leading competing qPCR Master Mixes. All amplifications were performed in accordance with the manufacturer's instructions. Hi-Tech Green 2X qPCR Universal Mix exhibited to be effective (mean Ct : 21.29) and highly specific because no second amplification signal could be identified with the melting curve. The results show that it has the best performance as compared with other qPCR Master Mix suppliers.

Materials Required but Not Provided

A compatible real-time PCR instrument, Vortex Mixer, Microcentrifuge, Plates and seals for your instruments

Procedure

Step 1: Prepare the PCR master mix

- Thaw Hi-Tech Green 2X qPCR Universal Master Mix and the rest of frozen reaction components to a temperature of 4°C. In order to entirely collect solutions, combine thoroughly and centrifuge briefly, then store at 4°C and avoid from light.
- Prepare (on ice or at room temperature) enough assay Master Mix for all reactions by adding all necessary components, except the DNA template, according to the recommendations in Table 1 (below).

Components	20 µl reaction	Final Conc.
PCR-grade water	Up to 20 µl	N/A
Hi-Tech Green 2X Mix	10 µl	1X
Forward Primer (10 µM)	0.8 µl or Variable	400 nM or Variable
Reverse Primer (10 µM)	0.8 µl or Variable	400 nM or Variable
Template DNA	As required	cDNA: 1pg–10ng Genomic DNA: 50ng-250ng or Variable

Step 2: Set up individual reactions

- Combine the assay Master Mix thoroughly to ensure consistency and equally dispense the solution into each qPCR tube or into the wells of a qPCR plate.
- Employ good pipetting practice to ensure assay precision and accuracy.
- Add DNA samples (and DNase-free H₂O if needed) to the PCR tubes or wells containing assay Master Mix (Table 1), seal the tubes or wells with flat caps or optically transparent film. Note: to ensure thorough mixing of reaction components, vortex for approximately 30 seconds (or more).
- Spin the tubes or plate to remove any air bubbles and collect the reaction mixture in the vessel bottom.
- Setup the thermal cycling protocol on a real-time PCR instrument according to below table. Note: optimization may be needed for better performance.
- Load the PCR tubes or plate into the real-time PCR instrument and commence the run.
- Perform data analysis according to the instrument-specific instructions.

Step 3: Run the PCR

Perform PCR with the following cycling protocol:

Step	Temperature	Duration	Cycle
Initial denaturation	95 °C	3-5 min ¹	1
Denaturation	95 °C	15 Sec	40
Annealing ² & Elongation	60 °C	60 Sec	
Melt analysis	OPTIONAL		

¹ Initial denaturation for 3 min at 95 °C is recommended for most assays. For GC-rich targets (>65% GC), 5 min at 95 °C may be used.

² An annealing temperature 5 °C lower than the calculated melting temperature (T_m) of the primer set is recommended as a first approach. If low yields and/or nonspecific amplification is obtained, an annealing temperature gradient PCR is recommended to determine the optimal annealing temperature of the primer pair.

Troubleshooting

Trouble	Cause	Solution
Poor Signal or No Signal	Inhibitor Present	Perform a dilution series of the PCR template to determine whether the effect of the inhibitory agent can be reduced. Take extra care with the nucleic acid extraction steps to minimize carryover of PCR inhibitors.
	Degraded Template Material	Do not store diluted template in water or at low concentrations. Check the integrity of template material by automated or manual gel electrophoresis.
	Inadequate Thermal Cycling Conditions	Try using a minimum extension time of 30 sec for genomic DNA and 15 sec for cDNA.
Signal in Negative Control	Contamination of Reaction Components with Target Sequence	To minimize the possibility of contamination of PCR components by PCR product or other template, designate a work area exclusively for PCR assay setup. Use a solution of 10% bleach instead of ethanol to prepare the workstation area for PCR assay setup. Ethanol will only induce precipitation of DNA in your work area, while the 10% bleach solution will hydrolyze, as well as dissolve, any residual DNA.
Poor Reproducibility Across Replicate Samples	Inhibitor Present	Perform a dilution series of the PCR template to determine whether the effect of the inhibitory agent can be reduced. Take extra care with the nucleic acid extraction steps to minimize carryover of PCR inhibitors.
	Primer Design	Verify primers design at different annealing temperatures.
Low or High Reaction Efficiency	Primer- Dimer	Reduce primer concentration. Evaluate primer sequences for complementarity and secondary structure. Redesign primers if necessary. Perform melt-curve analysis to determine if primer- dimers are present.
	Insufficient Optimization	Use a thermal gradient to identify the optimal thermal cycling conditions for a specific primer set.