

Glutathione (GSH) Reduced Assay Kit

Cat No. GSH-M-50

Size: 50 Reactions

Store at 4°C under dark conditions

MOLEQULE-ON®

Components	Quantity
Standard	10 mg
Reagent I	50 ml
Reagent II	50 ml
Reagent III	15 ml

Description

Glutathione (GSH) is a natural three peptide composed of glutamic acid (Glu), cysteine (Cys) and glycine (Gly), which is a compound containing sulfhydryl group (SH), it is an important antioxidant in plants, animals, fungi, and some bacteria and archaea. Glutathione reaction with 5,5'-dithiobis-2-nitrobenoic acid (DTNB) to form a yellow product. It has the maximum absorption at 412nm.

Reagents and Equipment Required but Not Provided:

Analytical balance, micro-homogenizer (2ml), refrigerated centrifuge, water-bath, Micropipette, Spectrophotometer, 1ml cuvette.

Sample Preparation

1. Tissue sample

Wash the fresh tissues twice with PBS, then add 0.1g washed tissue sample into homogenizer (Before using the homogenizer, wash with Solution I pre-cooled on ice). Add 1ml Solution I (the proportion of tissue and reagents can be kept constant), fully grind the tissue on ice (grinding under liquid nitrogen will have a better grinding effect). Centrifuge at 8000 rpm for 10 min at 4°C. Keep supernatant at 4 °C for test (The supernatant can be stored at -80°C for 10 days.)

2. Blood sample

a. Plasma: Centrifuge the sample at 600 g for 10 min at 4°C. Transfer the upper plasma into another tube and add same volume of Solution I, centrifuge at 8000 g for 10 min at 4°C. Place supernatant at 4°C for test. (The Supernatant can be stored at -80°C for 10 days).

b. Blood cell: Centrifuge the sample at 600 g for 10 min at 4°C. Discard the upper plasma, wash the settled blood cells with 3 times volumes of PBS (mix blood cell with PBS, centrifuge at 600 g for 10 min). Repeat the washing step two more times. In the washed blood cells, add equal volume of Solution I. Centrifuge at 8000 g for 10 min. Keep supernatant at 4°C for test (The supernatant can be stored at -80°C for 10 days).

3. Cell sample

Harvested cell should not be less than 10⁸. Wash with PBS twice (mix cell with PBS, centrifuge at 600 g for 10 min), wash with 3 times volumes of PBS for 3 times. Perform repeated freezing and thawing 2-3 times (suggest freezing in liquid nitrogen, and thawing at 37 °C water bath). Centrifuge at 8000 g for 10min. Place supernatant at 4 °C for test (The supernatant can be stored at -80°C for 10 days).

Protocol

1. Prewarm Spectrophotometer for 30 min, adjust the wavelength to 412 nm, set the counter to zero with distilled water.
2. Preheat Solution II in water bath: 37°C (mammal cell) , 25°C (other species).
3. Blank tube detection: Add 100 µl distilled water, 700 µl Solution II and 200 µl Solution III into 1 ml cuvette. Fully mixed for 2 min then determine absorbance A1 at 412nm.
4. Standard curve: 1 mg Standard dissolved with 1 ml distilled water (1mg/ml). Dilute Solution I tenfold, take appropriate solution to prepare standard solution with concentration of 200µg/ml, 100µg/ml, 50µg/ml, 25µg/ml, 12.5µg/ml. Add 100 µl Standard Solution, 700 µl Solution II , 200 µl Solution III into 1.5 ml centrifuge tube. Fully mixed and place for 2 min, detect the absorbance at 412nm, the values obtained subtract A1 as abscissa, according to the absorbance (x) and concentration (y, g/ml) to make the standard curve.
5. Sample determination: Add 100 µl sample, 700 µl Solution II , 200 µl Solution III into 1 ml cuvette. Fully mixed and place for 2 min, detect absorbance as A2 at 412nm, $\Delta A = A2 - A1$.

Calculations

According to the standard curve, sample ΔA into the formula (x), calculate the sample concentration of Y (g/ml).

1. Protein concentration

$$\text{GSH}(\mu\text{g} / \text{mg protein}) = y * V_{rv} / V_s / C_{pr} = y * 10 / C_{pr}$$

2. Sample weight

$$\text{GSH}(\mu\text{g} / \text{g}) = y * V_{rv} / (V_s / V_{sv} * W) = y * 10 / W$$

3. Cell amount

$$\text{GSH}(\mu\text{g} / 10^4 \text{ cell}) = y * V_{rv} / (V_s / V_{sv} * n) = y * 10 / n$$

4. Solution volume

$$\text{GSH}(\mu\text{g} / \text{mL}) = y * V_{rv} / V_s = y * 10$$

n: Cell amount

V_{rv} : Total reaction volume, 1ml;

V_{sv} : Total supernatant volume, 1 ml;

V_s : The volume of supernatant was added into the reaction system, 100µl=0.1 ml;

W: Sample weight, g;

C_{pr} : Supernatant protein concentration, mg/ml.

Notes

1. The sample needs to be treated completely. If the test cannot be completed temporarily, it can be stored at -80°C.
2. The standard solution should be prepared just before use.
3. If the content of GSH in the sample is uncertain, dilute the sample to several gradients for test.
4. Solution I contained protein precipitant, the supernatant could not be used for protein concentration determination.