

Microcentrifuge Tubes

CAT No.	DESCRIPTION
MQ-BS-0.5-1000	0.5 ml Microcentrifuge tube, natural, bag of 1000
MQ-BS-1.5-1000	1.5 ml Microcentrifuge tube, natural, bag of 1000
MQ-BS-2.0-1000	2.0 ml Microcentrifuge tube, natural, bag of 1000

MOLEQULE-ON®



Description

MOLEQULE-ON is committed to provide high quality, low retention centrifuge tubes. These 1.5ml and 2.0ml microcentrifuge tubes are manufactured from FDA approved resin and free from RNase, DNase, DNA and Pyrogen. Leak proof seal lids to facilitate single hand operation. Frosted area is available for labeling. No sharp points or edges, ultra smooth design that does not catch gloves. Clear and precise graduation with immaculate visibility to perform all types of experiments. These are come in natural color and can be centrifuged at up to 20,000 RCF. The tubes are safe to use at temperatures ranges from - 80 to +121°C. Can be fit in all standard and high capacity centrifuge rotors. These are sterile thus suitable for any routine and sensitive laboratory and diagnostic work. 1.5ml tubes come in bulk pack of 1000 while 2.0ml size is available in 500 tubes packing. In addition, 1.5 ml tubes attributed to take total 1.7 ml liquid, this feature minimize the chances of 1.5 ml liquid to expel out when the tube lid will be closed.