


Large Electrophoresis Tub

Cat No: HGE-M-003

MOLEQULE-ON[®]



Description

The Large Tub from MOLEQULE-ON provides a wide range of versatility. Three different size gel trays are available 15 X 7cm, 15 X 10cm and 15 x 15cm allowing the selection of one, two or all three gel length according to the experiment. Also, additional purchases of accessories are no longer required. This wide range of tray options allow up to 210 samples to be resolved per gel. The large size gel of 15cm run length allows restriction fragment or other close molecular weight sample bands to be easily separated and identified. Speed loading is accomplished using 10, 14, 16, 18, 28 and 30 sample multi-channel pipette compatible combs.

Features

- Three gel tray options.
- Run up to 210 samples.
- Low buffer volumes.
- Multichannel pipette compatible combs for speed loading.



All 3 tray sizes for optimum flexibility and value


Specifications

Unit Dimensions (W x L x H)	17.5 x 26.5 x 9cm
Gel Dimensions (W x L)	15 X 7cm
	15 X 10cm
	15 X 15cm
Sample Capacity	15 x 7cm Tray- 70 Samples
	15 X 10cm Tray - 140 Samples
	15 X 15cm Tray - 210 Samples
Buffer Volume	500ml

Additional Accessories (on request)

- Combs with colour code & sample number 1, 2, 4, 10, 10MC*, 12, 14MC, 16, 18MC, 20, 28MC, 30MC, 35
- Thicknesses 0.75, 1, 1.5, 2MM


* Multi Channel Pipette Compatible Combs



Buffer saver blocks Physically reduce the volume of a gel chamber and so reduce buffer requirements, saving cost.



Scoops allow safe transfer of gels



Loading guides allow easy well identification transfer of gel and sample loading



Combs with different colour codes for different thicknesses

Accessories

Cat No:	Description
MQ-MS15-UV7	Large Tub, 15 x 7cm UV Tray
MQ-MS15-UV10	Large Tub, 15 x 10cm UV Tray
MQ-MS15-UV15	Large Tub, 15 x 15cm UV Tray
MQ-MS15-UVDAM	Large Tub, UV Tray Dams
MQ-MS15-PE	Large Tub - Positive Electrode
MS15-NE	Large Tub - Negative Electrode
MQ-MS15-LG	Loading Guides
MQ-MS15-WP	Viewing Platform
MQ-MSCHOICECP	Cool-pack and platform for Large Tub
MQ-MSCHOICEBSB	Buffer Saver Blocks for Large Tub pack/2 saves 190ml of buffer
MQ-MS15/20-FC	Large Tub Flexi caster
MQ-MS15-UVS	15cm UV Gel Scoop
MQ-MS15-1-0.75	Comb Prep 1, Marker 1, 0.75mm thick

MQ-MS15-2-0.75	Comb Prep 2, Marker 2, 0.75mm thick
MQ-MS15-4-0.75	Comb Prep 4, Marker 2, 0.75mm thick
MQ-MS15-10-0.75	Comb 10 sample, 0.75mm thick
MQ-MS15-12-0.75	Comb 12 sample, 0.75mm thick
MQ-MS15-20-0.75	Comb 20 sample, 0.75mm thick
MQ-MS15-35-0.75	Comb 35 sample, 0.75mm thick
MQ-MS15-10MC-0.75	Comb 10 sample MC, 0.75mm thick
MQ-MS15-14MC-0.75	Comb 14 sample MC, 0.75mm thick
MQ-MS15-16MC-0.75	Comb 16 sample MC, 0.75mm thick
MQ-MS15-18MC-0.75	Comb 18 sample MC, 0.75mm thick
MQ-MS15-28MC-0.75	Comb 28 sample MC, 0.75mm thick
MQ-MS15-30MC-0.75	Comb 30 sample MC, 0.75mm thick
MQ-MS15-1-1	Comb Prep 1, Marker 1, 1mm thick
MQ-MS15-2-1	Comb Prep 2, Marker 2, 1mm thick
MS15-4-1	Comb Prep 4, Marker 2, 1mm thick
MQ-MS15-10-1	Comb 10 sample, 1mm thick
MQ-MS15-12-1	Comb 12 sample, 1mm thick
MQ-MS15-20-1	Comb 20 sample, 1mm thick
MQ-MS15-35-1	Comb 35 sample, 1mm thick
MQ-MS15-10MC-1	Comb 10 sample MC, 1mm thick
MQ-MS15-14MC-1	Comb 14 sample MC, 1mm thick
MQ-MS15-16MC-1	Comb 16 sample MC, 1mm thick
MQ-MS15-18MC-1	Comb 18 sample MC, 1mm thick
MQ-MS15-28MC-1	Comb 28 sample MC, 1mm thick
MQ-MS15-30MC-1	Comb 30 sample MC, 1mm thick
MQ-MS15-1-1.5	Comb Prep 1, Marker 1, 1.5mm thick
MQ-MS15-2-1.5	Comb Prep 2, Marker 2, 1.5mm thick
MQ-MS15-4-1.5	Comb Prep 4, Marker 2, 1.5mm thick
MQ-MS15-10-1.5	Comb 10 sample, 1.5mm thick
MQ-MS15-12-1.5	Comb 12 sample, 1.5mm thick
MQ-MS15-20-1.5	Comb 20 sample, 1.5mm thick
MQ-MS15-35-1.5	Comb 35 sample, 1.5mm thick
MQ-MS15-30MC-1.5	Comb 30 sample MC, 1.5mm thick
MQ-MS15-10MC-1.5	Comb 10 sample MC, 1.5mm thick
MQ-MS15-14MC-1.5	Comb 14 sample MC, 1.5mm thick
MQ-MS15-16MC-1.5	Comb 16 sample MC, 1.5mm thick
MQ-MS15-18MC-1.5	Comb 18 sample MC, 1.5mm thick
MQ-MS15-28MC-1.5	Comb 28 sample MC, 1.5mm thick
MQ-MS15-1-2	Comb Prep 1, Marker 1, 2mm thick
MQ-MS15-2-2	Comb Prep 2, Marker 2, 2mm thick
MQ-MS15-4-2	Comb Prep 4, Marker 2, 2mm thick
MQ-MS15-10-2	Comb 10 sample, 2mm thick
MQ-MS15-12-2	Comb 12 sample, 2mm thick
MQ-MS15-16MC-2	Comb 16 sample MC, 2mm thick
MQ-MS15-20-2	Comb 20 sample, 2mm thick
MQ-MS15-30MC-2	Comb 30 sample MC, 2mm thick
MQ-MS15-10MC-2	Comb 10 sample MC, 2mm thick
MQ-MS15-14MC-2	Comb 14 sample MC, 2mm thick
MQ-MS15-18MC-2	Comb 18 sample MC, 2mm thick
MQ-MS15-28MC-2	Comb 28 sample MC, 2mm thick
MQ-MS15-35-2	Comb 35 sample, 2mm thick
MS15-x -yLG	Combination comb and Loading guide for choice horizontal, X = comb well number, y = comb thickness

Additional Accessories (on request)

- Combs with colour code & sample number 1, 2, 4, 8MC*, 8, 10, 12MC and 16
- Thicknesses 0.75, 1, 1.5, 2MM

* Multi Channel Pipette Compatible Combs



Buffer saver blocks Physically reduce the volume of a gel chamber and so reduce buffer requirements, saving cost.



Scoops allow safe transfer of gels



Loading guides allow easy well identification transfer of gel and sample loading



Combs with different colour codes for different thicknesses