

Micropipettes (Fixed Volume)



Cat No.	Range
MPF-M-5	5 µl
MPF-M-10	10 µl
MPF-M-20	20 µl
MPF-M-25	25 µl
MPF-M-50	50 µl
MPF-M-100	100 µl
MPF-M-200	200 µl
MPF-M-250	250 µl
MPF-M-500	500 µl
MPF-M-1000	1000 µl
MPF-M-2000	2000 µl
MPF-M-5000	5000 µl
MPF-M-10000	10000 µl



MOLEQULE-ON[®]

Description

MOLEQULE-ON pipette sets a new standard for safety and ergonomics in mechanical pipettes. It minimizes the risk of Repetitive Stress Injury (RSI) by requiring lower force for tip attachment and detachment. The micropipette is fully autoclavable. Wide range of micropipette sizes starting from 5µl up to 10 mL is available. This ergonomic pipette has user-friendly handle design with soft grip for comfortable operation to prevent hand stress ensure ultimate comfort.

Features

- 13 volume micropipette models with range from 5 µl to 10 ml.
- UV resistant.
- Ergonomic design with soft grip
- Fully autoclavable for 15 minutes at 121°C.
- Suitable for most brands of pipette tips.
- Volume ranges colour coded.
- Volumes adjustable, with easy-to-read display.