

Description

The display is fully visible during pipette operation with good size digits, and the large soft handle pad ensures easy tip ejection. The micropipette is fully autoclavable. Wide range of micropipette sizes starting from 0.1μ l up to 10 mL is available. This ergonomic pipette has user- friendly handle design with soft grip for comfortable operation to prevent hand stress ensure ultimate comfort.

Features

- > Nine variable volume micropipette models with range from 0.1 μ l to 10 ml.
- ➢ UV resistant.
- Ergonomic design with soft grip
- ➤ Fully autoclavable for 15 minutes at 121°C.
- Suitable for most brands of pipette tips.
- Volume ranges colour coded.
- Volumes adjustable, with easy-to-read display.