

# Multichannel Micropipettes



**MOLEQULE-ON®**

8 Channel		12 Channel	
Range	Cat No.	Range	Cat No.
0.5-10ul	MCP-M-8-10	0.5-10ul	MCP-M-12-10
2-20ul	MCP-M-8-20	2-20ul	MCP-M-12-20
5-50ul	MCP-M-8-50	5-50ul	MCP-M-12-50
10-100ul	MCP-M-8-100	10-100ul	MCP-M-12-100
20-200ul	MCP-M-8-200	20-200ul	MCP-M-12-200
30-300ul	MCP-M-8-300	30-300ul	MCP-M-12-300

## Description

MOLEQULE-ON offers Multichannel Micropipettes, which make pipetting easy and effortless. MQ Multichannel pipettes are completely autoclavable without dismantling it. This Multichannel Pipette is made with superior quality material for better protection against chemical and physical corrosion. It has enhanced UV resistance. The very light weight model makes it easy to operate even for extended time. Fatigue free working with convenient finger rest. Tip cone is designed for smooth fitting with many brands of tips. In house cleaning, repair and calibration is possible without any special tools.

## Features

- Innovative Spring and seal design for one of the lowest plunging force requirements
- All internal component in the manifold is made of high-quality alloys produced by fine tooling to deliver precise volume across all channel
- Easily removable manifold for convenient maintenance
- 360° rotatable manifold for effortless and flexible twisting
- Highly durable PVDF tip cone
- Different size with different color code
- Universal tip compatible mean will accept wide range of tips
- Available in 8 & 12 channel



For lesser strain & proper tip fitment for precise pipetting



For consistently precise results