

Bottle Top Dispenser

Cat No: BDR-M-10

MOLEQULE-ON®



Description

MOLEQULE-ON's Bottle Top Dispenser (1-10ml) offers unmatched features and satisfaction to customer. It is made from absolute best quality raw material for precision and durability that last for many years. It is developed with a unique springless valve design that offers high-chemical resistance and longer life with minimal maintenance.

Features

- Telescoping filling tube, compatible with a bottle-height range of 170mm to 300mm
- Discharge tube with 360 rotation for convenient dispensing in multiple bottles
- Single tool recalibration
- Ergonomic design for convenient sampling and low pressure dispensing for large volumes
- Snap-fit mounting of discharge tube for easy handling
- Safety locking feature for prevention against accidental spillage
- Multiple sizes are available

