

Constant Temperature and Humidity Incubator

Cat No: THI-M-160

MOLEQULE-ON[®]



Description

MQ Constant Temperature and Humidity Incubator provides a temperature range from 0°C to 65°C and humidity ranges from 50-90% RH. The 160L chamber volume (comes in various sizes under different catalogue number) offers ample space to meet customer requirement. Chamber can be sterilized by built-in UV lamp for new startup.

Features

- The starting up time can be preset.
- Universal castors for easy movement.
- Microprocessor temperature control.
- CFC-free refrigerant with high efficiency.
- Double door design and magnetic door sealing.
- UV lamp for sterilization.
- The working chamber adopts round angle structure, easy to clean.
- Over temperature protection, temperature sensor malfunction alarm.
- High and low water level alarm, automatic heating off for low water level.
- Parameters can be saved automatically when power off or system error.

Specifications

Model	SIN-M-12000
Capacity	160L
Temp. Range	0~65°C
Temp. Fluctuation	±1 °C
Temp. Stability/Uniformity	±2 °C
Humidity Range	50~90% RH (±5~8% RH)
Sterilization Type	UV lamp
Shelves	3 pcs/adjustable
Consumption	710W
Power Supply	220V, 50/60Hz
Internal Size (W*D*H)mm	450*420*850
External Size (W*D*H)mm	590*620*1400