

MQ Shaking Incubator (Vertical Type)

Cat No: SIN-M-12000

MOLEQULE-ON[®]



Description

MQ Shaking Incubator (Vertical Type) is a large incubator which can accommodate number of flasks and bottles to facilitate shaking and incubation of different samples. It offers PID controller with LCD display. Forced convection circulation with rotary shaking mode. It can work at temperature of 5°C to 60°C continuously up to 500 hour.

Features

- Audio and visual alarm for over temperature.
- Alarm for over-speed
- Leakage or over-current protection

Specifications

Model	SIN-M-12000
Controller	PID
Display	LCD
Circulation Mode	Forced convection
Shaking Mode	Rotary shaking
Drive Mode	Multidimensional drive
Rotating Speed	40~300rpm
Rotating Accuracy	±1rpm
Temp. Range	RT+5~60°C
Temp. Accuracy	±0.1°C
Temp. Uniformity	±1°C(At 37°C)
Timing Range	0~500 hours
Shaking Plate Size	500*350mm
Shaking Plate Quantity	2pcs
Max. Capacity	1000ml*12 or 500ml*24 or 250ml*30 or 100ml*48 or 50ml*40