

Mini Centrifuge

Cat No: MCM-M-6000

MOLEQULE-ON[®]



Description

MOLEQULE-ON introduces new easy-to-change SnapSpin™ rotors in its mini centrifuge. Two new SnapSpin rotors (a six-place standard rotor and an eight-place strip rotor) Snap-in and out-for fast and easy change without using any tool. Two sets of additional adapters accommodate all major micro tubes and PCR tube styles (2.0, 1.5, 0.5, 0.2 mL) for customized use.

This mini centrifuge is ideal for micro filter cell separations, molecular biology applications and HPLC samples. With most spins lasting around 20 seconds, the ergonomically designed lid is specifically designed to fit in the palm of the hand for easy repeat opening and use. Load samples and for repeated spins, simply close the lid to reactivate.

Features

- Faster operation with Snap-in rotor.
- Machine will not operate with the lid open.
- Made with UV resistant plastic
- Small footprint offers minimal bench-top space and removable cord for compact storage.

MOLEQULE-ON Mini Centrifuge includes:

- Two easy changeable, tool-free rotors:
 - Six-place standard rotor for 1.5/2 ml tubes
 - Eight-place strip rotor for two 0.2ml tube strips or 16 single tubes
- Six 0.2 ml and six 0.5 ml adapters for customized use
- Storage Case – Keeps rotors and adapters contained and dust-free
- Tube Rack – Convenient for handling samples before and after quick spin



Tube Racks

Storage Case

Rotors and Adapters

Specifications

Maximum speed:	6,000 RPM
Maximum RCF	2000 x g
Dimensions	153 x 128 x 104mm
Warranty	One-year warranty