


# MQ Multi-well Plates

Ordering Info		
Cat No.	Product	Packaging
CCP-M-6	6 Well cell culture plate, PS, TC Treated, Sterile	Box of 100 Plates
CCP-M-12	12 Well cell culture plate, PS, TC Treated, Sterile	Box of 100 Plates
CCP-M-24	24 Well cell culture plate, PS, TC Treated, Sterile	Box of 100 Plates
CCP-M-48	48 Well cell culture plate, PS, TC Treated, Sterile	Box of 100 Plates
CCP-M-96	96 Well cell culture plate, PS, TC Treated, Flat Bottom, Sterile	Box of 100 Plates
	Inserts for 6, 12, 24 and 48 well plates	



MOLEQULE-ON offers different multi-well plates to meet the demand of every laboratory. These Plates are delicately design to prevent cross-contamination. These are free from RNase/Dnase and pyrogen. The air-vent covers facilitate gas exchange and homogeneous cell culture. Alphanumeric identification for easy samples locating. 96 wells plate applicable to HTS (High Throughput Screening).