

PCR Plates

Cat No.	Product	Packaging
PP96-NS-0.1-N	MQ PCR Plate, 96 well, Non-Skirted, 0.1ml, Natural	15 pcs/box, Case of 10 boxes
PP96-NS-0.1-W	MQ PCR Plate, 96 well, Non-Skirted, 0.1ml, White	15 pcs/box, Case of 10 boxes
PP96-NS-0.1-B	MQ PCR Plate, 96 well, Non-Skirted, 0.1ml, Black	15 pcs/box, Case of 10 boxes
PP96-SS-0.1-N	MQ PCR Plate, 96 well, Semi-Skirted, 0.1ml, Natural	15 pcs/box, Case of 10 boxes
PP96-SS-0.1-W	MQ PCR Plate, 96 well, Semi-Skirted, 0.1ml, White	15 pcs/box, Case of 10 boxes
PP96-NS-0.2-N	MQ PCR Plate, 96 well, Non-Skirted, 0.2ml, Natural	15 pcs/box, Case of 10 boxes
PP96-NS-0.2-W	MQ PCR Plate, 96 well, Non-Skirted, 0.2ml, White	15 pcs/box, Case of 10 boxes
PP96-NS-0.2-B	MQ PCR Plate, 96 well, Non-Skirted, 0.2ml, Black	15 pcs/box, Case of 10 boxes
PP96-SS-0.2-N	MQ PCR Plate, 96 well, Semi-Skirted, 0.2ml, Natural	15 pcs/box, Case of 10 boxes
PP96-SS-0.2-W	MQ PCR Plate, 96 well, Semi-Skirted, 0.2ml, White	15 pcs/box, Case of 10 boxes
PP96-SS-0.2-B	MQ PCR Plate, 96 well, Semi-Skirted, 0.2ml, Black	15 pcs/box, Case of 10 boxes
PP96-FS-0.2-N	MQ PCR Plate, 96 well, Full-Skirted, 0.2ml, Natural	15 pcs/box, Case of 10 boxes

MOLEQULE-ON®



Description

MQ PCR Plates are DNase/RNase free and non-pyrogenic. Plates come in two sizes with Non-Skirted, Semi-Skirted and Full-skirted format. They are available in three colors to meet the demand of every lab. The innovative product is designed for general and quantitative PCR applications, suitable for most regular PCR machines such as ABI/LT, AGILENT, BIORAD, EPPENDORF etc. The plates can be cut into 16-well, 24-well, 32-well or 48-well plates for using together with 8-strip caps. The sealing films ensures greater excellency during thermal transfer.