


# Serological Pipettes



The serological pipettes by MOLEQULE-ON are sterile, individually wrapped and non-pyrogenic. They have both descending & ascending black graduations marks. Top grade resins reduces liquid binding and improves sampling accuracy. Filter plug prevents the liquid from entering pipetting device and their capacity calibration rate is  $\pm 2\%$ . They are available in six volumes.

Ordering Info		
Cat No.	Product	Packaging
SP-M-1	1 ml serological pipette	50/unit, 20units/case
SP-M-2	2 ml serological pipette	50/unit, 20units/case
SP-M-5	5 ml serological pipette	50/unit, 4units/case
SP-M-10	10 ml serological pipette	50/unit, 20units/case
SP-M-25	25 ml serological pipette	25/unit, 8units/case
SP-M-50	50 ml serological pipette	25/unit, 4units/case