

# Technical Data Sheet

# MOLEQULE-ON<sup>®</sup>

## Fetal Equine Serum (Horse Serum)

Cat #: TC-M-FES01200

### **Description:**

Horse Serum can be used as an alternative to fetal bovine serum for growth factors and hormones in cell culture studies. It is mostly used as a blocking agent in immunoassays. Horse Serum contains more immunoglobulins than fetal bovine serum (FBS). It has an increased protein content, and is a cost-effective alternative to fetal bovine serum.

### **Features:**

- Alternative to Fetal Bovine Serum
- Suitable for ES cell culture
- Autologous cell culture system

### **Storage:**

Temperature: -20°C

*The product is for general laboratory use. Not for human diagnostic or therapeutic use.*

# MOLEQULE-ON