

MEM Non-Essential Amino Acids Solution TC-M-NEAA100 100X Size: 100ml

Description

Contains the same non-essential amino acids found in the standard Minimum Essential Medium (MEM) at a strength of 100X. It is used as a growth supplement for cell culture medium, to increase cell growth and viability.

Formulation

Components	Molecular Weight	Concentration (mg/L)	mM
Amino Acids			
Glycine	75	750	10
L-Alanine	89	890	10
L-Asparagine	132	1320	10
L-Aspartic acid	133	1330	10
L-Glutamic Acid	147	1470	10
L-Proline	115	1150	10
L-Serine	105	1050	10

Storage

Store at 2-8°C.

For laboratory use only. Not for human diagnostic and therapeutic use.