

## Sterile Water, Nuclease Free Cat No. CH-M-120



## **DESCRIPTION**

Nuclease free water offered by MOLEQULE-ON, New Zealand is made to use in all kinds of molecular biology applications. It has been prepared without the use of chemicals such as DEPC (diethylpyrocarbonate) using an in-house method. It is membrane-filtered and tested for DNase and RNase activity.

## **STORAGE**

Store at 18-25°C.

For laboratory use only. Not intended for any animal or human therapeutic or diagnostic uses.

P.O. Box: 15549, New Lynn, Auckland: 0640, New Zealand Phone: +64 27 3759195 E-mail: info@molegule-on.com, Web: www.molegule-on.com