

Triton X-100 Detergent Solution

Cat No. BF-M-TXT-10

Description

Triton X-100 is a commonly used detergent in various laboratories. Triton X-100 Detergent Solution contain 10% highly purified Triton X-100 which is used as non-ionic detergent for the isolation of membrane proteins under non-denaturing conditions. It is used a multi-purpose detergent to optimized protein-protein interactions.

STORAGE

Store at 15-25°C.

For laboratory use only. Not intended for human therapeutic or diagnostic use.

MOLEQULE-ON