

# Technical Data Sheet

**MOLEQULE-ON<sup>®</sup>**

## Tris Buffer, 1M (pH 8.0)

Cat No. BF-M-TB-1M

### Description

Tris (hydroxymethyl)-aminomethane hydrochloride buffer in 1M concentration is the most suitable solution to be used in a variety of biological applications. The buffer is mostly used to control pH in different molecular biology procedures.

### STORAGE

Store at 15-25°C.

For laboratory use only. Not intended for human therapeutic or diagnostic use.

MOLEQULE-ON