

# Technical Data Sheet

**MOLEQULE-ON<sup>®</sup>**

## **Phosphate Buffer, 1M**

**Cat No. BF-M-PB-1M**

### **Description**

Phosphate buffers (also called Gomori buffers) consist of Potassium phosphate monobasic and Potassium phosphate dibasic salts. The buffer is soluble in water and has high buffering capacity. Phosphate buffer is used in different research domain including molecular biology.

### **STORAGE**

Store at 15-25°C.

For laboratory use only. Not intended for human therapeutic or diagnostic use.

MOLEQULE-ON