


Cell Culture Dishes and Flasks



MOLEQULE-ON cell culture dishes for microscopic applications are made of optically clear polymers. They have enhanced lid design that ensures easy gas exchange. External grip enables easy grasping. These dishes are stackable. Sterile, RNase/DNase free and non-pyrogenic.

Ordering Info		
Cat No.	Product	Packaging
CCD-M-35	Cell Culture Dishes, 35mm, TC Treated, PS, Sterile	20 pcs/pack, 25 packs/case
CCD-M-60	Cell Culture Dishes, 60mm, TC Treated, PS, Sterile	20 pcs/pack, 25 packs/case
CCD-M-100	Cell Culture Dishes, 100mm, TC Treated, PS, Sterile	10 pcs/pack, 20 packs/case



MQ Cell culture flasks are ergonomically designed for superior handling while achieving lowest level of contamination. These have short flat bevel angle which enable easy access. They can be stackable for excellent space utilization. All the cell culture consumables are DNase/RNase free and Non-pyrogenic. All flasks models are available with plugged and filtered caps.

Ordering Info		
Cat No.	Product	Packaging
TF-M-TC-25	25cm ² cell culture flask; PS, 40ml, breathable cap, TC treatment, sterilization	5pcs/pack, 40 bags/carton
TF-M-TC-75	75cm ² cell culture flask; PS, 170ml, breathable cap, TC treatment, sterilization	5pcs/pack, 20 bags/carton
TF-M-TC-175	175cm ² cell culture flask; PS, 450ml, breathable cap, TC treatment, sterilization	5pcs/pack, 10 packs/carton