

# MQ Ultrasonic Cleaner

Cat No: USC-M-30

**MOLEQULE-ON**<sup>®</sup>



## Description

MQ Ultrasonic Cleaner is a multifunctional device with large LCD display to show the precise time, temperature and power level. It includes a drain valve and humanized carrying handle. A build-in “sound insulation paste” keeps the machine with low noise. It comes with sweep and degas function. The sleeping mode is available to keep electrical energy saving.

## Features

- Ultra-low noise, adjustable time, adjustable power, frequency sweep function, degassing function, memory function, sleeping mode.
- Time setting with 1-99 minutes and built-in heating to 60°C , adjustable power range of 10%~100%.
- Memory function will automatically save the most recently used setup parameters for next use.

## Specifications

<b>Model</b>	<b>USC-M-30</b>
Ultrasonic Power (W) Max	600
Ultrasonic Power Setting (%)	10%-100% (step by 10%)
Heating Power (W) Max	500
Capacity (L) Max	30
External Size (LxWxH) cm	53x33x39
Transducer Qty (PCS)	10
Ultrasonic Frequency (KHz)	40
Temperature setting (°C) Max	Default 60°C
Timer Setting (Min)	1-99