

Magnetic Beads and Retriever Bar



Magnetic beads offered by MOLEQULE-ON are cylindrical stirring bars encapsulated with high-grade, crack-resistant PTFE (polytetrafluoroethylene) coating. They are used in stirring small to large volume of liquid or viscous mixtures. They have a smooth surface and central pivot ring that ensure efficient spinning even on curved or uneven bases. The stirring bars are resistant to all chemical and temperatures up to 260°C. The retriever bar can easily take out the magnetic bead from any solution.



Ordering Info		
Product	Size	Cat Number
Magnetic Beads	Small	MBS-M-5
Magnetic Beads	Medium	MBM-M-5
Magnetic Beads	Large	MBL-M-5
Retriever Bar		RMB-M-01